

Anti-ENPP5 Monoclonal Antibody-Biotin

Product Details

| | |
|-------------------|-----------------------------|
| Ig Type: | Rabbit monoclonal IgG |
| Reactivity: | Human |
| Conjugation: | Biotin |
| Molecular Weight: | 150 kDa |
| Purification: | Protein A Affinity Purified |

Applications

| | |
|--------------------|---|
| Verified Activity: | Flow cytometry analysis of ENPP5 overexpressed 465F cells with TMAZ-0271B, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line). |
| Application: | ELISA, FCM |
| Recommended | 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA |

Properties

| | |
|----------------------|---|
| Purity: | > 95% as determined by SDS-PAGE. |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|--|
| Immunogen: | ENPP5 |
| Antigen Species: | Human |
| Gene ID: | 59084 |
| Uniprot ID: | Q9UJA9 |
| Synonyms: | Ectonucleotide pyrophosphatase/phosphodiesterase family member 5 |
| Biology Area: | Other Functions |

Research Background

Can hydrolyze NAD but cannot hydrolyze nucleotide di- and triphosphates. Lacks lysopholipase D activity. May play a role in neuronal cell communication

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481

Anti-ENPP5 Monoclonal Antibody-Biotin

Product Details

Molecular Weight: 150 kDa

Properties

Antigen Details

Synonyms: Ectonucleotide pyrophosphatase/phosphodiesterase family member 5

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481