

Anti-SHISAL1 Monoclonal Antibody-Biotin

Product Details

Molecular Weight: 150 kDa

Properties

Antigen Details

Synonyms: Protein shisa-like-1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481

Anti-SHISAL1 Monoclonal Antibody-Biotin

Product Details

| | |
|-------------------|-----------------------------|
| Ig Type: | Rabbit monoclonal IgG |
| Reactivity: | Human |
| Conjugation: | Biotin |
| Molecular Weight: | 150 kDa |
| Purification: | Protein A Affinity Purified |

Applications

| | |
|--------------------|---|
| Verified Activity: | Flow cytometry analysis of SHISAL1 overexpressed 484F cells with TMAZ-0276B, followed by goat anti-rabbit IgG-ABflo 647 (red line). The isotype control is rabbit IgG (black line). |
| Application: | ELISA,FCM |
| Recommended | 0.1-0.2 µg/10E6 cells for FCM; 1 ng/µl for ELISA |

Properties

| | |
|----------------------|---|
| Purity: | > 95% as determined by SDS-PAGE. |
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|----------------------|
| Immunogen: | SHISAL1 |
| Antigen Species: | Human |
| Gene ID: | 85352 |
| Uniprot ID: | Q3SXP7 |
| Synonyms: | Protein shisa-like-1 |
| Biology Area: | Other Functions |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481