

## Compound T129533

## Chemical Properties


CAS No. :

Formula: C<sub>19</sub>H<sub>22</sub>O<sub>7</sub>

Molecular Weight: 362.38

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

4-ethyl-7-hydroperoxy-10-methoxy-3,5,12-trimethyl-14-oxatetracyclo[9.2.1.0<sup>1,5</sup>.0<sup>7,11</sup>]tetradeca-3,9-diene-2,6,8-trione is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7595 mL	13.7977 mL	27.5953 mL
5 mM	0.5519 mL	2.7595 mL	5.5191 mL
10 mM	0.276 mL	1.3798 mL	2.7595 mL
50 mM	0.0552 mL	0.276 mL	0.5519 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481