

HIV-IN-3

## Chemical Properties

CAS No. :

Formula: C<sub>21</sub>H<sub>32</sub>ClN<sub>7</sub>O<sub>3</sub>

Molecular Weight: 465.98

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	HIV-IN-3 (Compound 22a) is a potent HIV inhibitor (IC <sub>50</sub> : 1.5 μM) with potential for studying HIV-related diseases.
Targets(IC <sub>50</sub> )	Others,HIV Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.146 mL	10.7301 mL	21.4601 mL
5 mM	0.4292 mL	2.146 mL	4.292 mL
10 mM	0.2146 mL	1.073 mL	2.146 mL
50 mM	0.0429 mL	0.2146 mL	0.4292 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481