


## Sulfadoxine-D4

## Chemical Properties

CAS No. :	1330266-05-1
Formula:	C <sub>12</sub> H <sub>10</sub> D <sub>4</sub> N <sub>4</sub> O <sub>4</sub> S
Molecular Weight:	314.35
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sulfadoxine-D4 is a deuterated compound of Sulfadoxine. Sulfadoxine (T0836) has a CAS number of 2447-57-6. Sulfadoxine (T0836) is a broad-spectrum sulfanilamide and a synthetic analog of para-aminobenzoic acid (PABA) with bacteriostatic and antimalarial properties.
Targets(IC50)	HIV Protease,Endogenous Metabolite,Antibiotic,Parasite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1812 mL	15.9058 mL	31.8117 mL
5 mM	0.6362 mL	3.1812 mL	6.3623 mL
10 mM	0.3181 mL	1.5906 mL	3.1812 mL
50 mM	0.0636 mL	0.3181 mL	0.6362 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Raman J, et al. Am J Trop Med Hyg. 2010, 82(5):788-94.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481