


DL-Epinephrine-D3

Chemical Properties

CAS No. :	1189977-29-4
Formula:	C ₉ H ₁₀ D ₃ NO ₃
Molecular Weight:	186.22
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	DL-Epinephrine-D3 is a deuterated compound of DL-Epinephrine. DL-Epinephrine (T1071L2) has a CAS number of 329-65-7. DL-Epinephrine (DL-Adrenaline) is a hormone and a neurotransmitter.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.370 mL	26.850 mL	53.6999 mL
5 mM	1.074 mL	5.370 mL	10.740 mL
10 mM	0.537 mL	2.685 mL	5.370 mL
50 mM	0.1074 mL	0.537 mL	1.074 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481