


Fluconazole-D4

Chemical Properties

CAS No. :	1124197-58-5
Formula:	C ₁₃ H ₈ D ₄ F ₂ N ₆ O
Molecular Weight:	310.3
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fluconazole-D4 is a deuterated compound of Fluconazole. Fluconazole (T1388) has a CAS number of 86386-73-4. Triazole antifungal agent that is used to treat oropharyngeal candidiasis and cryptococcal meningitis in AIDS.
Targets(IC50)	Antibiotic, Antifungal

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2227 mL	16.1134 mL	32.2269 mL
5 mM	0.6445 mL	3.2227 mL	6.4454 mL
10 mM	0.3223 mL	1.6113 mL	3.2227 mL
50 mM	0.0645 mL	0.3223 mL	0.6445 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481