


Azathioprine-D3

Chemical Properties

CAS No. :	2702733-53-5
Formula:	C9D3H4N7O2S
Molecular Weight:	280.28
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azathioprine-D3 is a deuterated compound of Azathioprine. Azathioprine (T1237) has a CAS number of 446-86-6. Azathioprine, an immunosuppressive agent, inhibits purine synthesis and GTP-binding protein Rac1 activation.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5679 mL	17.8393 mL	35.6786 mL
5 mM	0.7136 mL	3.5679 mL	7.1357 mL
10 mM	0.3568 mL	1.7839 mL	3.5679 mL
50 mM	0.0714 mL	0.3568 mL	0.7136 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481