

Tetrabromobisphenol A-D4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Tetrabromobisphenol A-D4 is a deuterated compound of Tetrabromobisphenol A. Tetrabromobisphenol A (T71753) has a CAS number of 79-94-7.
-------------	-----------------------------------------------------------------------------------------------------------------------------------------

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481