


Dehydroepiandrosterone-2,2,3,4,4,6-D6

Chemical Properties

CAS No. :	1261254-39-0
Formula:	C ₁₉ H ₂₂ D ₆ O ₂
Molecular Weight:	294.46
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dehydroepiandrosterone-2,2,3,4,4,6-D6 is a deuterated compound of Dehydroepiandrosterone. Dehydroepiandrosterone (T0856L) has a CAS number of 53-43-0. Dehydroepiandrosterone (T0856L) (DHEA) is an endogenous steroid hormone that is secreted primarily by the adrenal gland and is the most abundant sex steroid. Although the mechanisms of action of exogenously administered DHEA have not been fully illuminated, they may result in both direct and indirect physiologic effects. Direct effects include GABA-a receptor complex and NMDA receptor modulation and enhanced pancreatic beta cell insulin secretion and anti-glucocorticoid activities.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.396 mL	16.9802 mL	33.9605 mL
5 mM	0.6792 mL	3.396 mL	6.7921 mL
10 mM	0.3396 mL	1.698 mL	3.396 mL
50 mM	0.0679 mL	0.3396 mL	0.6792 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481