


Tetrahydro 11-Deoxycorticosterone-D3

Chemical Properties

CAS No. : 72205-58-4
 Formula: C₂₁H₃₁D₃O₃
 Molecular Weight: 337.51
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tetrahydro 11-Deoxycorticosterone-D3 is a deuterated compound of Tetrahydro 11-Deoxycorticosterone. Tetrahydro 11-Deoxycorticosterone (T37725) has a CAS number of 567-03-3.
Targets(IC50)	GABA Receptor, Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9629 mL	14.8144 mL	29.6288 mL
5 mM	0.5926 mL	2.9629 mL	5.9258 mL
10 mM	0.2963 mL	1.4814 mL	2.9629 mL
50 mM	0.0593 mL	0.2963 mL	0.5926 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481