

## Nandrolone 17-propionate-D3

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Nandrolone 17-propionate-D3 is a deuterated compound of Nandrolone 17-propionate. Nandrolone 17-propionate (T19470) has a CAS number of 7207-92-3. Nandrolone propionate is an anabolic steroid action primarily consisting of androgenic and nitrogen-retentive properties.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481