

## 20(S)-Hydroxy Prednisolone-D4

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

20(S)-Hydroxy Prednisolone-D4 is a deuterated compound of 20 (S)-Hydroxy Prednisolone. 20 (S)-Hydroxy Prednisolone has a CAS number of 2299-46-9.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481