


## Levofloxacin-D8

### Chemical Properties

CAS No. : 1217716-71-6  
 Formula: C<sub>18</sub>H<sub>12</sub>D<sub>8</sub>FN<sub>3</sub>O<sub>4</sub>  
 Molecular Weight: 369.42  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	Levofloxacin-D8 is a deuterated compound of Levofloxacin. Levofloxacin (T6567) has a CAS number of 100986-85-4. Levofloxacin (T6567) is a broad-spectrum antibiotic topoisomerase II and topoisomerase IV inhibitor, used to treat respiratory, urinary tract, gastrointestinal, and abdominal infections.
Targets(IC50)	Antibacterial, Antibiotic

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7069 mL	13.5347 mL	27.0695 mL
5 mM	0.5414 mL	2.7069 mL	5.4139 mL
10 mM	0.2707 mL	1.3535 mL	2.7069 mL
50 mM	0.0541 mL	0.2707 mL	0.5414 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481