


Salbutamol-D9

Chemical Properties

CAS No. :	1173021-73-2
Formula:	C ₁₃ H ₁₂ D ₉ NO ₃
Molecular Weight:	248.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Salbutamol-D9 is a deuterated compound of Salbutamol. Salbutamol (T1139) has a CAS number of 18559-94-9. Albuterol stimulates beta2-adrenergic receptors in the lungs, thereby activating the enzyme adenylate cyclase that catalyzes the conversion of ATP to cyclic-3', 5'-adenosine monophosphate (cAMP). Albuterol Sulfate is the sulfate salt of the short-acting sympathomimetic agent albuterol, a 1:1 racemic mixture of (R)-albuterol and (S)-albuterol with bronchodilator activity. Increased cAMP concentrations relax the bronchial smooth muscle, relieve bronchospasms, and reduce inflammatory cell mediator release, especially from mast cells. To a lesser extent, albuterol stimulates beta1-adrenergic receptors, thereby increasing the force and rate of myocardial contraction.
Targets(IC50)	Adrenergic Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0263 mL	20.1313 mL	40.2625 mL
5 mM	0.8053 mL	4.0263 mL	8.0525 mL
10 mM	0.4026 mL	2.0131 mL	4.0263 mL
50 mM	0.0805 mL	0.4026 mL	0.8053 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481