


Lacosamide-D3

Chemical Properties

CAS No. : 1217689-95-6
 Formula: C₁₃H₁₅D₃N₂O₃
 Molecular Weight: 253.31
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lacosamide-D3 is a deuterated compound of Lacosamide. Lacosamide (T2540) has a CAS number of 175481-36-4. Lacosamide is a functionalized amino acids. Lacosamide modulates sodium channels in a novel manner: It selectively enhances sodium channel slow inactivation with no effects on fast inactivation Lacosamide has demonstrated antiepileptic effectiveness in different rodent seizure models and antinociceptive potential in experimental animal models that reflect distinct types and symptoms of neuropathic as well as chronic inflammatory pain.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9477 mL	19.7387 mL	39.4773 mL
5 mM	0.7895 mL	3.9477 mL	7.8955 mL
10 mM	0.3948 mL	1.9739 mL	3.9477 mL
50 mM	0.079 mL	0.3948 mL	0.7895 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481