


Stavudine- $\alpha,\alpha,\alpha,6$ -D4

Chemical Properties

CAS No. :	1219803-67-4
Formula:	C10H8D4N2O4
Molecular Weight:	228.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Stavudine- $\alpha,\alpha,\alpha,6$ -D4 is a deuterated compound of Stavudine. Stavudine (T1404) has a CAS number of 3056-17-5. Stavudine, a nucleoside reverse transcriptase inhibitor analog of thymidine, has activity against HIV.
Targets(IC50)	Apoptosis,Nucleoside Antimetabolite/Analog,HIV Protease,NOD-like Receptor (NLR), Reverse Transcriptase,Autophagy

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3814 mL	21.9068 mL	43.8135 mL
5 mM	0.8763 mL	4.3814 mL	8.7627 mL
10 mM	0.4381 mL	2.1907 mL	4.3814 mL
50 mM	0.0876 mL	0.4381 mL	0.8763 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481