

Cinoxacino-D8

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Cinoxacino-D8 is a deuterated compound of Cinoxacino. Cinoxacino (T0274) has a CAS number of 28657-80-9. Cinoxacin, an older synthetic antimicrobial, was related to the quinolone class of antibiotics. Its activity was similar to oxolinic acid and nalidixic acid.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481