


Keto Pioglitazone-D4

Chemical Properties

CAS No. :	1215370-26-5
Formula:	C ₁₉ H ₁₄ D ₄ N ₂ O ₄ S
Molecular Weight:	374.45
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Keto Pioglitazone-D4 is a deuterated compound of Keto Pioglitazone. Keto Pioglitazone (T20465) has a CAS number of 146062-45-5. Ketopioglitazone is an active metabolite of pioglitazone.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6706 mL	13.3529 mL	26.7058 mL
5 mM	0.5341 mL	2.6706 mL	5.3412 mL
10 mM	0.2671 mL	1.3353 mL	2.6706 mL
50 mM	0.0534 mL	0.2671 mL	0.5341 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481