


## Adipic Acid-D10

## Chemical Properties

CAS No. :	25373-21-1
Formula:	C <sub>6</sub> H <sub>10</sub> O <sub>4</sub>
Molecular Weight:	156.203
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Adipic Acid-D10 is a deuterated compound of Adipic Acid. Adipic Acid (T2228) has a CAS number of 124-04-9. Adipic acid is a precursor for the production of nylon.
Targets(IC50)	Endogenous Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.402 mL	32.0102 mL	64.0205 mL
5 mM	1.2804 mL	6.402 mL	12.8041 mL
10 mM	0.6402 mL	3.201 mL	6.402 mL
50 mM	0.128 mL	0.6402 mL	1.2804 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481