


Atorvastatin-D5 Calcium Salt

Chemical Properties

CAS No. :	222412-82-0
Formula:	C ₆₆ H ₅₈ D ₁₀ CaF ₂ N ₄ O ₁₀
Molecular Weight:	1165.40
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Atorvastatin-D5 Calcium Salt is a deuterated compound of Atorvastatin Calcium Salt. Atorvastatin Calcium Salt (T3116) has a CAS number of 134523-00-5.
Targets(IC50)	HMG-CoA Reductase, Autophagy

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8581 mL	4.2904 mL	8.5807 mL
5 mM	0.1716 mL	0.8581 mL	1.7161 mL
10 mM	0.0858 mL	0.429 mL	0.8581 mL
50 mM	0.0172 mL	0.0858 mL	0.1716 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481