

## rac-Aprepitant-D4 (enantiomer)

## Chemical Properties


CAS No. :

Formula: C<sub>23</sub>H<sub>17</sub>D<sub>4</sub>F<sub>7</sub>N<sub>4</sub>O<sub>3</sub>

Molecular Weight: 538.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	rac-Aprepitant-D4 (enantiomer) is a deuterated compound of rac-Aprepitant. rac-Aprepitant has a CAS number of 170729-80-3.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8572 mL	9.2859 mL	18.5718 mL
5 mM	0.3714 mL	1.8572 mL	3.7144 mL
10 mM	0.1857 mL	0.9286 mL	1.8572 mL
50 mM	0.0371 mL	0.1857 mL	0.3714 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481