

Melitracen-D6 hydrochloride

Chemical Properties


CAS No. :

Formula: C₂₁H₂₀D₆ClN

Molecular Weight: 333.93

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Melitracen-D6 hydrochloride is a deuterated compound of Melitracen Hydrochloride (TMSM-1551).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9946 mL	14.9732 mL	29.9464 mL
5 mM	0.5989 mL	2.9946 mL	5.9893 mL
10 mM	0.2995 mL	1.4973 mL	2.9946 mL
50 mM	0.0599 mL	0.2995 mL	0.5989 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481