

Bromocriptine-13C-D3

Chemical Properties


CAS No. :

Formula: 13CC31D3H37BrN5O5

Molecular Weight: 658.61

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

Bromocriptine-13C-D3 the 13C and deuterated compound of Bromocriptine. Bromocriptine has a CAS number of 25614-03-3. Bromocriptine is a semisynthetic ergotamine alkaloid. It is a dopamine D2 agonist. It suppresses prolactin secretion.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5183 mL	7.5917 mL	15.1835 mL
5 mM	0.3037 mL	1.5183 mL	3.0367 mL
10 mM	0.1518 mL	0.7592 mL	1.5183 mL
50 mM	0.0304 mL	0.1518 mL	0.3037 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481