

Esomeprazole-D3 Sodium Salt

Chemical Properties


CAS No. :

Formula: C17H15D3N3NaO3S

Molecular Weight: 370.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Esomeprazole-D3 Sodium Salt is a deuterated compound of Esomeprazole Sodium Salt (T2686L).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6996 mL	13.4982 mL	26.9964 mL
5 mM	0.5399 mL	2.6996 mL	5.3993 mL
10 mM	0.270 mL	1.3498 mL	2.6996 mL
50 mM	0.054 mL	0.270 mL	0.5399 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481