


Nitrobenzene-D5

Chemical Properties

CAS No. :	4165-60-0
Formula:	C6D5HONO2
Molecular Weight:	128.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Nitrobenzene-D5 is a deuterated compound of Nitrobenzene.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.804 mL	39.0198 mL	78.0396 mL
5 mM	1.5608 mL	7.804 mL	15.6079 mL
10 mM	0.7804 mL	3.902 mL	7.804 mL
50 mM	0.1561 mL	0.7804 mL	1.5608 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481