


1,4-Dichlorobenzene-D4

Chemical Properties

CAS No. : 3855-82-1
 Formula: C₆D₄H₀Cl₂
 Molecular Weight: 151.03
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,4-Dichlorobenzene-D4 is a deuterated compound of 1,4-Dichlorobenzene. 1,4-Dichlorobenzene (T7431) has a CAS number of 106-46-7. 1,4-Dichlorobenzene (T7431) is a disinfectant, pesticide, and deodorant.
-------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6212 mL	33.106 mL	66.212 mL
5 mM	1.3242 mL	6.6212 mL	13.2424 mL
10 mM	0.6621 mL	3.3106 mL	6.6212 mL
50 mM	0.1324 mL	0.6621 mL	1.3242 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481