


rac 1-Palmitoyl-3-Chloropropanediol-D5

Chemical Properties

CAS No. : 1346599-60-7
 Formula: C₁₉H₃₂D₅ClO₃
 Molecular Weight: 353.98
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	rac 1-Palmitoyl-3-Chloropropanediol-D5 is a deuterated compound of rac 1-Palmitoyl-3-Chloropropanediol. rac 1-Palmitoyl-3-Chloropropanediol (TMSM-4319) has a CAS number of 30557-04-1.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.825 mL	14.1251 mL	28.2502 mL
5 mM	0.565 mL	2.825 mL	5.650 mL
10 mM	0.2825 mL	1.4125 mL	2.825 mL
50 mM	0.0565 mL	0.2825 mL	0.565 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481