


Biphenyl-D10

Chemical Properties

CAS No. :	1486-01-7
Formula:	C12D10
Molecular Weight:	164.27
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Biphenyl-D10 is a deuterated compound of Biphenyl (PDK0398). Biphenyl (PDK0398) has a CAS number of 92-52-4.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.0875 mL	30.4377 mL	60.8754 mL
5 mM	1.2175 mL	6.0875 mL	12.1751 mL
10 mM	0.6088 mL	3.0438 mL	6.0875 mL
50 mM	0.1218 mL	0.6088 mL	1.2175 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481