


## Triphenyl-D15 Phosphate

## Chemical Properties

CAS No. : 1173020-30-8  
 Formula: C18D15O4P  
 Molecular Weight: 341.38  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Triphenyl-D15 Phosphate is a deuterated compound of Triphenyl Phosphate. Triphenyl Phosphate (T203211) has a CAS number of 115-86-6.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9293 mL	14.6464 mL	29.2929 mL
5 mM	0.5859 mL	2.9293 mL	5.8586 mL
10 mM	0.2929 mL	1.4646 mL	2.9293 mL
50 mM	0.0586 mL	0.2929 mL	0.5859 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481