


2, 2'-Dipyridyl-D8

Chemical Properties

CAS No. : 32190-42-4
 Formula: C10D8N2
 Molecular Weight: 164.23
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2, 2'-Dipyridyl-D8 is a deuterated compound of 2, 2'-Dipyridyl. 2, 2'-Dipyridyl (Fr14395) has a CAS number of 366-18-7.
-------------	-------------------------------------------------------------------------------------------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.089 mL	30.4451 mL	60.8902 mL
5 mM	1.2178 mL	6.089 mL	12.178 mL
10 mM	0.6089 mL	3.0445 mL	6.089 mL
50 mM	0.1218 mL	0.6089 mL	1.2178 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481