


Sulfameter-D4

Chemical Properties

CAS No. :	1189483-96-2
Formula:	C11H8D4N4O3S
Molecular Weight:	284.33
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Sulfameter-D4 is a deuterated compound of Sulfameter. Sulfameter (T0748) has a CAS number of 651-06-9. Sulfameter (T0748) is a long-acting sulfonamide antibacterial agent used to treat urinary tract infections and leprosy.
Targets(IC50)	Antibacterial, Antibiotic

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.517 mL	17.5852 mL	35.1704 mL
5 mM	0.7034 mL	3.517 mL	7.0341 mL
10 mM	0.3517 mL	1.7585 mL	3.517 mL
50 mM	0.0703 mL	0.3517 mL	0.7034 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481