


## N-Methyltryptamine-D3

## Chemical Properties

CAS No. : 1794756-39-0  
 Formula: C<sub>11</sub>H<sub>11</sub>D<sub>3</sub>N<sub>2</sub>  
 Molecular Weight: 177.26  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-Methyltryptamine-D3 is a deuterated compound of N-Methyltryptamine. N-Methyltryptamine has a CAS number of 61-49-4.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6414 mL	28.2072 mL	56.4143 mL
5 mM	1.1283 mL	5.6414 mL	11.2829 mL
10 mM	0.5641 mL	2.8207 mL	5.6414 mL
50 mM	0.1128 mL	0.5641 mL	1.1283 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481