


Pyrene-D10

Chemical Properties

CAS No. :	1718-52-1
Formula:	C16D10
Molecular Weight:	212.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pyrene-D10 is a deuterated compound of Pyrene. Pyrene (T200798) has a CAS number of 129-00-0.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.7101 mL	23.5505 mL	47.1009 mL
5 mM	0.942 mL	4.7101 mL	9.4202 mL
10 mM	0.471 mL	2.355 mL	4.7101 mL
50 mM	0.0942 mL	0.471 mL	0.942 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481