

CAM4066

## Chemical Properties


CAS No. : 2101206-81-7

Formula: C<sub>26</sub>H<sub>26</sub>ClN<sub>3</sub>O<sub>4</sub>

Molecular Weight: 479.96

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CAM4066 is a novel potent and selective CK2alpha inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0835 mL	10.4175 mL	20.8351 mL
5 mM	0.4167 mL	2.0835 mL	4.167 mL
10 mM	0.2084 mL	1.0418 mL	2.0835 mL
50 mM	0.0417 mL	0.2084 mL	0.4167 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481