

Perforatumone

Chemical Properties

CAS No. :	1600490-96-7
Formula:	C35H52O5
Molecular Weight:	552.79
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Perforatumone is a useful organic compound for research related to life sciences and the catalog number is T124148.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.809 mL	9.045 mL	18.0901 mL
5 mM	0.3618 mL	1.809 mL	3.618 mL
10 mM	0.1809 mL	0.9045 mL	1.809 mL
50 mM	0.0362 mL	0.1809 mL	0.3618 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481