

2-(dec-9-en-1-yl)-3-hydroxypentanedioic acid

Chemical Properties

CAS No. :

Formula: C₁₅H₂₆O₅

Molecular Weight: 286.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	2-(dec-9-en-1-yl)-3-hydroxypentanedioic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.492 mL	17.4599 mL	34.9199 mL
5 mM	0.6984 mL	3.492 mL	6.984 mL
10 mM	0.3492 mL	1.746 mL	3.492 mL
50 mM	0.0698 mL	0.3492 mL	0.6984 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481