

Compound T130493

Chemical Properties


CAS No. :

Formula: C₃₄H₄₂O₁₉

Molecular Weight: 754.69

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-(((2R,3S,4S,5R,6S)-3,4-dihydroxy-6-((5-hydroxy-2-(4-hydroxy-3-methoxyphenyl)-4-oxo-3,4-dihydro-2H-1-benzopyran-7-yl)oxy)-5-(((2S,3R,4R,5R,6S)-3,4,5-trihydroxy-6-methyloxan-2-yl)oxy)oxan-2-yl)methoxy)-3-hydroxy-3-methyl-5-oxopentanoic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.325 mL	6.625 mL	13.2505 mL
5 mM	0.265 mL	1.325 mL	2.6501 mL
10 mM	0.1325 mL	0.6625 mL	1.325 mL
50 mM	0.0265 mL	0.1325 mL	0.265 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481