

## [(2R,3S,4S,5R,6S)-3,4,5-trihydroxy-6-(3-hydroxy-5-methylphenoxy)oxan-2-yl]methyl trihydroxybenzoate

## Chemical Properties


CAS No. :

Formula: C<sub>20</sub>H<sub>22</sub>O<sub>11</sub>

Molecular Weight: 438.39

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	[(2R,3S,4S,5R,6S)-3,4,5-trihydroxy-6-(3-hydroxy-5-methylphenoxy)oxan-2-yl]methyl 3,4,5-trihydroxybenzoate is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2811 mL	11.4054 mL	22.8107 mL
5 mM	0.4562 mL	2.2811 mL	4.5621 mL
10 mM	0.2281 mL	1.1405 mL	2.2811 mL
50 mM	0.0456 mL	0.2281 mL	0.4562 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481