

DSTYLSSTLTLSK acetate

Chemical Properties

CAS No. :


Formula: C66H111N15O28

Molecular Weight: 1562.67

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	DSTYLSSTLTLSK acetate is a human generic peptide that can be used for the quantitative detection of infliximab.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.6399 mL	3.1997 mL	6.3993 mL
5 mM	0.128 mL	0.6399 mL	1.2799 mL
10 mM	0.064 mL	0.320 mL	0.6399 mL
50 mM	0.0128 mL	0.064 mL	0.128 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Jing-Ya Shiao, et al. Development of efficient on-bead protein elution process coupled to ultra-high performance liquid chromatography-tandem mass spectrometry to determine immunoglobulin G subclass and glycosylation for discovery of bio-signatures in pancreatic disease. *Chromatogr A*. 2020 Jun 21;1621:461039.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481