

Super-TDU

Chemical Properties

CAS No. : 1599441-71-0

Formula: C237H370N66O69S

Molecular Weight: 5279.94

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

SVDDHFAKSLGDTWLQIGSGNPKTANVPQTVPMRLRKLPSDFKPPPE

Biological Description

Description	Super-TDU, an inhibitory peptide, targets YAP-TEADs interaction.
Targets(IC50)	Others,YAP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.1894 mL	0.947 mL	1.894 mL
5 mM	0.0379 mL	0.1894 mL	0.3788 mL
10 mM	0.0189 mL	0.0947 mL	0.1894 mL
50 mM	0.0038 mL	0.0189 mL	0.0379 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481