


2,16-Kauranediol 2-O-beta-D-allopyranoside

Chemical Properties

CAS No. : 474893-07-7
 Formula: C₂₆H₄₄O₇
 Molecular Weight: 468.62
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2,16-Kauranediol 2-O-beta-D-allopyranoside is a natural product for research related to life sciences. The catalog number is TN2675 and the CAS number is 474893-07-7.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1339 mL	10.6696 mL	21.3393 mL
5 mM	0.4268 mL	2.1339 mL	4.2679 mL
10 mM	0.2134 mL	1.067 mL	2.1339 mL
50 mM	0.0427 mL	0.2134 mL	0.4268 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481