


Hexaethylene glycol dimethyl ether

Chemical Properties

CAS No. : 1072-40-8

Formula: C14H30O7

Molecular Weight: 310.387



Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Hexaethylene glycol dimethyl ether, a PEG-based PROTAC linker, facilitates the synthesis of PROTACs.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, composed of two different ligands connected by a linker—one targeting an E3 ubiquitin ligase and the other the target protein—utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2218 mL	16.1088 mL	32.2175 mL
5 mM	0.6444 mL	3.2218 mL	6.4435 mL
10 mM	0.3222 mL	1.6109 mL	3.2218 mL
50 mM	0.0644 mL	0.3222 mL	0.6444 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481