

RKI-1447

Chemical Properties


CAS No. : 1342278-01-6

Formula: C₁₆H₁₄N₄O₂S

Molecular Weight: 326.37

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	RKI-1447 is a potent inhibitor of ROCK1 and ROCK2. It has anti-invasive and antitumor activities.
Targets(IC50)	ROCK
In vitro	RKI-1447 inhibits the growth of 87% of breast tumors, reducing the average volume of breast tumors by 7.7 times compared to the control mice. In a transgenic mouse model, RKI-1447 effectively suppresses breast tumor growth. When tumors in mice are treated with RKI-1447, the average percentage increase in tumor volume is only 8.8%.
In vivo	RKI-1447 inhibits the phosphorylation of ROCK substrates MLC-2 and MYPT-1 in human cancer cells, without affecting the phosphorylation levels of AKT, MEK, and S6 kinases at concentrations up to 10 μM. It suppresses the migration, invasion, and anchorage-independent growth of breast cancer cells. The crystal structure of the RKI-1447/ROCK1 complex reveals that RKI-1447 is a type I kinase inhibitor, binding to the ATP site through interactions with the hinge region and DFG motif.
Kinase Assay	Z-Lyte FRET kinase assay: Kinase inhibition is measured using the Invitrogen Z-Lyte? FRET kinase assay with Ser/Thr 13 peptide substrate based on the myosin light chain sequence KKRPPQRRYSNVF. Compounds are tested on three separate days with 8 point dilutions performed in duplicate to determine average IC50 values. The assay conditions are optimized to 15 μL of kinase reaction volume with 5 ng of enzyme in 50 mM HEPES (pH 7.5), 10 mM MgCl ₂ , 1 mM EGTA, and 0.01% Brij-35. The reaction is incubated for 1 h at room temperature in the presence of 1.5 μM of peptide substrate with 12.5 μM of ATP (for ROCK1) or 2 μM of substrate with 50 μM of ATP (for ROCK2). The reaction is then stopped and the ratio of phosphorylated to unphosphorylated peptides is determined by selective cleavage of only the unphosphorylated peptide as described by the manufacturer. This is followed by excitation of coumarin at 400 nm resulting in emission at 445 nm and energy transfer to fluorescein and final emission at 520 nm. The substrate contains both coumarin and fluorescein and only uncleaved phosphorylated substrate will undergo FRET. The ratio of the signals at 445 nm and 520 nm is measured using a Wallac EnVision Plate Reader, model 2102 plate-reader.
Cell Research	Cells are plated in a 96 well tissue culture plate (1200 cells per well) and incubated for 24 hours. After incubation the cells are treated with vehicle or increasing concentrations of RKI-1447 for 72 hours. After incubation, freshly prepared MTT (2 mg/ml) is added to each well and incubated for 3 hours. After incubation the plates are read at 540 nm.

A DRUG SCREENING EXPERT

Cell Research	(Only for Reference)
---------------	----------------------

Solubility Information

Solubility	DMSO: 18.33 mg/mL (56.16 mM),Sonication is recommended. H2O: < 1 mg/mL (insoluble or slightly soluble), Ethanol: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (6.13 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.064 mL	15.320 mL	30.6401 mL
5 mM	0.6128 mL	3.064 mL	6.128 mL
10 mM	0.3064 mL	1.532 mL	3.064 mL
50 mM	0.0613 mL	0.3064 mL	0.6128 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Patel RA, et al. Cancer Res, 2012, 72(19), 5025-5034.

Pireddu R, et al. Medchemcomm, 2012, 3(6), 699-709.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481