


## Securisteroside

## Chemical Properties

CAS No. :	54964-57-7
Formula:	C <sub>51</sub> H <sub>88</sub> O <sub>7</sub>
Molecular Weight:	813.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Securisteroside is a natural product for research related to life sciences. The catalog number is TN5324 and the CAS number is 54964-57-7.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2296 mL	6.1482 mL	12.2965 mL
5 mM	0.2459 mL	1.2296 mL	2.4593 mL
10 mM	0.123 mL	0.6148 mL	1.2296 mL
50 mM	0.0246 mL	0.123 mL	0.2459 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481