


Phenglutarimid

Chemical Properties

CAS No. :	1156-05-4
Formula:	C17H24N2O2
Molecular Weight:	288.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Phenglutarimid is an anticholinergic. It is used as an antiparkinsonian agent.
Targets(IC50)	Others,AChR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4676 mL	17.3382 mL	34.6765 mL
5 mM	0.6935 mL	3.4676 mL	6.9353 mL
10 mM	0.3468 mL	1.7338 mL	3.4676 mL
50 mM	0.0694 mL	0.3468 mL	0.6935 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

BATTEGAY R. [Aturban (phenglutarimide), a new therapeutic drug for parkinsonism; comparative investigation on the therapy of parkinsonism].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481