


GW-475151

## Chemical Properties

CAS No. : 221323-43-9  
 Formula: C<sub>19</sub>H<sub>28</sub>N<sub>4</sub>O<sub>5</sub>  
 Molecular Weight: 424.51  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	GW-475151 inhibits human neutrophil elastase (HNE).
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3557 mL	11.7783 mL	23.5566 mL
5 mM	0.4711 mL	2.3557 mL	4.7113 mL
10 mM	0.2356 mL	1.1778 mL	2.3557 mL
50 mM	0.0471 mL	0.2356 mL	0.4711 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481