


Viscidulin III tetraacetate

Chemical Properties

CAS No. :	96684-81-0
Formula:	C ₂₅ H ₂₂ O ₁₂
Molecular Weight:	514.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Viscidulin III tetraacetate is a natural product for research related to life sciences. The catalog number is TN5239 and the CAS number is 96684-81-0.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9439 mL	9.7193 mL	19.4386 mL
5 mM	0.3888 mL	1.9439 mL	3.8877 mL
10 mM	0.1944 mL	0.9719 mL	1.9439 mL
50 mM	0.0389 mL	0.1944 mL	0.3888 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481