


Orilotimod potassium

Chemical Properties

CAS No. :	960155-19-5
Formula:	C ₁₆ H ₁₈ KN ₃ O ₅
Molecular Weight:	371.43
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Orilotimod potassium (APO805K1) can be used to study plaque psoriasis.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 60 mg/mL (161.54 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6923 mL	13.4615 mL	26.923 mL
5 mM	0.5385 mL	2.6923 mL	5.3846 mL
10 mM	0.2692 mL	1.3461 mL	2.6923 mL
50 mM	0.0538 mL	0.2692 mL	0.5385 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Tam TF, et al. Processes for the manufacture of pharmaceutically acceptable salts of thymodepressin. WO2008074128 A1 .

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481