


(+)-Hannokinol

Chemical Properties

CAS No. :	408324-76-5
Formula:	C ₁₉ H ₂₄ O ₄
Molecular Weight:	316.39
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(+)-Hannokinol is a natural product for research related to life sciences. The catalog number is TN6383 and the CAS number is 408324-76-5.
Targets(IC50)	NO Synthase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1607 mL	15.8033 mL	31.6066 mL
5 mM	0.6321 mL	3.1607 mL	6.3213 mL
10 mM	0.3161 mL	1.5803 mL	3.1607 mL
50 mM	0.0632 mL	0.3161 mL	0.6321 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Shim K S, Hwang Y H, Jang S A, et al. Ethanol Extract of Amomum tsao-ko Ameliorates Ovariectomy-Induced Trabecular Loss and Fat Accumulation[J]. Molecules. 2021, 26(4): 784.

Shim K S, Hwang Y H, Jang S A, et al. Ethanol Extract of Amomum tsao-ko Ameliorates Ovariectomy-Induced Trabecular Loss and Fat Accumulation. Molecules. 2021 Feb 3;26(4):784. doi: 10.3390/molecules26040784.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481